

SYSTEM APPLICATION




TPG[®] AND

TPG[®] uno

IMMEDIATE LOADING
DENTAL IMPLANT SYSTEM

SIMPLADENT[®]









TWO-PIECE IMPLANTS

TPG®



TPG® implants provide a sharp and cutting crestal thread, which engages in to the corticals of the extraction sockets. The apical compression thread provides perfect stability both in compressed spongy and in cortical bone. A large variety of abutments is available. **TPG®** implants may be used both in extraction sockets and in healed bone areas. They are used for circular bridges, segments, and single tooth replacement. Under adequate loading conditions and if good primary stability is given, the treatment may be performed in an immediate load protocol. **TPG®** implants are not indicated for single tooth replacements, except if two or more implants provide antirotation for the implant-to-prosthetic construction when replacing one tooth. **TPG®** implants may be used in combination with other compression screw implants (e.g. KOC®) as a base for fixed prosthetics.

TPG® IMPLANTS II POLISHED, MATERIAL Ti6AL4V*

	Description	Max. Ø	Endosseous Ø	Endosseous length	Height above bone	REF	Price cat.
	TPG 4.1 8	4.95 mm	4.8 mm	8 mm	1.8 mm	13-420106	K
	TPG 4.1 10	4.95 mm	4.8 mm	10 mm	1.8 mm	13-420107	K
	TPG 4.1 12	4.95 mm	4.8 mm	12 mm	1.8 mm	13-420104	K
	TPG 4.1 15	4.95 mm	4.8 mm	15 mm	1.8 mm	13-420100	K
	TPG 4.1 17	4.95 mm	4.8 mm	17 mm	1.8 mm	13-420101	K
	TPG 4.1 19	4.95 mm	4.8 mm	19 mm	1.8 mm	13-420102	K
	TPG 4.1 21	4.95 mm	4.8 mm	21 mm	1.8 mm	13-420103	K
	TPG 4.1 23	4.95 mm	4.8 mm	23 mm	1.8 mm	13-420105	K

Internal Cone Connection
8 degr.

All **TPG®** implants are delivered incl. OF5 STI abutment (REF 13-420426)


TPG 4.1 8

TPG 4.1 10

TPG 4.1 12

TPG 4.1 15

TPG 4.1 17

TPG 4.1 19


TPG 4.1 21

TPG 4.1 23



***Ti6AL4VELI**, also called „Grade 5“, is a purified version of the conventional 6/4 Ti-Alloy. It is used for more than 50% of all human metallic implants. The material features superb mechanical and chemical properties, and it is very resistant against corrosion. Also regarding bio-compatibility and the support of cell-growth this specific titanium-alloy shows advantages compared to pure titanium. For these reasons the material is used so widespread. Ti6AL4V ELI provides a 25% higher stability compared to pure titanium.

EXAMPLE FOR THE USE OF TPG® IMPLANTS


TPG® implants may be inserted into fresh extraction sockets with lateral and bi-cortical engagement.








View from top: **TPG®** implants may be inserted into fresh extraction. The voids in the extraction sockets may receive augmentation or fill with woven bone.




TPG® implant engaged in the distal maxilla and the pterygoid plate of the sphenoid bone. In this indication this implant may be used for circular bridges and segments with at least 3 implants.

Please note: Additional abutments, tools and accessories can be found in the Simpladent® **Two-Part** system application (REF 13-0004-05)

ABUTMENTS FOR CEMENTED PROSTHETICS II WITH CONUS, WITHOUT ANTI-ROTATION






	Description	Code	REF	Price cat.
	Straight, screw-in abutment with pre-assembled screw for cemented crowns or bridges. Tighten with TT 1.25 .	TLA STI	13-425110	D
	15° angled, screw-in abutment with pre-assembled screw for cemented crowns or bridges. Tighten with TT 1.25 .	TLA 15 STI	13-425111	E
	20° angled, screw-in abutment with pre-assembled screw for cemented crowns or bridges. Tighten with TT 1.25 .	TLA 20 STI	13-425112	E
	25° angled, screw-in abutment with pre-assembled screw for cemented crowns or bridges. Tighten with TT 1.25 .	TLA 25 STI	13-425113	E
	One-piece abutment with solid thread, for cemented crowns or bridges, height above implant 6.5 mm, 8° incline. Can be shortened or trimmed. Flat on one side. Direct impression or transfer to lab analog IA STI. Tighten with HT 1.25 .	TCA STI	13-420313	D

UCLA ABUTMENTS

	Description	Code	REF	Price cat.
	Castable abutment incl. SF 215 replacement screw, 10.5 mm high. For ULCA. Tighten with HT 1.25 .	PA STI	13-420358	C

Use with impression post TSX STI (REF13-420411) and analog IA STI (REF13-420347)

PLATFORM OFFSET ABUTMENTS

	Description	Neck	Code	REF	Price cat.
	Platform offset abutment. Tool: IT K (REF 13-462320)*	1.5 mm	OF 4 STI	13-420424*	E
	Platform offset abutment, large head.	1.3 mm	OF 5 STI	13-420426	E
	Platform offset abutment, large head.	4 mm	OF 5L STI	13-420427	E
	Insertion tool, length 20 mm. Tool: RAT, RAT2, TW2		IT K	13-462320	C
	Insertion tool for large head		IT2 BCS	13-900030	E

* Laboratory accessories can be found in the Simpladent **Single-Piece** system application (REF 13-0003-05)




HEATLESS® DRILLS „DOS“ FOR TPG® IMPLANTS

Surgical steel, color coded, depth coded and autoclaveable. The drill is marked with laser depth markings.









	Description	max. working length	Code	REF	Price cat.
	BCD 1, yellow	15 mm	BCD 1	13-900240	C
	BCDX 1, yellow	15 mm	BCDX 1	13-900243	C
	DOS 1, yellow	17 mm	DOS 1	455311	D
	DOS 2, black	17 mm	DOS 2	455312	D
	DOS 3, red	17 mm	DOS 3	455313	D
	Twist Drill 2.0 / 21, Ø 2.0 mm	21 mm	Twist Drill 2.0	13-90022	D
	Countersink head, for implant-Ø4.1 mm		DRA 4.1 CS	13-425047	D





CORTICAL MILLING

	Description		Code	REF	Price cat.
	C-Drill KM1 3.0 - 3.2	Cortical milling	C-Drill KM1	13-455300	E
	C-Drill KM2 3.7 - 4.1	Cortical milling	C-Drill KM2	13-455301	E
	C-Drill KM3 5.0	Cortical milling	C-Drill KM3	13-455302	E











INSERTION TOOLS FOR RATCHET

	Description		Code	REF	Price cat.
	Insertion tool, short 12 mm		IT1 TPG	13-425116	K
	Insertion tool, medium 16 mm		IT2 TPG	13-425117	K
	Insertion tool, long 20 mm		IT3 TPG	13-425118	K
	Insertion tool for contra-angle, 30 mm		ITM TPG	13-425119	B
	Ratchet, for all Hex tools and insertion tools		RAT 2	13-425051	K
	Torque wrench 10 - 70 Ncm. It is recommended to have the torque ratchets recalibrated by us once a year.		TW2	13-425402	S

TOOLS FOR MANUAL INSERTION

	Description	Code	REF	Price cat.
	Engages in the internal trilobe inside the conus. Must be used with max. 60 Ncm only. If the bone is harder, IT1-3 TPG has to be used with RAT 2 or TW2.	Adapter TPG	13-425106	D
For two part TPG*				
	Adapter short, for large head	AA5 15°	13-462052	C
	Adapter short, for large head	AA5 25°	13-462053	C
	Handgrip, not demountable. Cleaning in ultrasonic bath at 45°, with alkaline cleansing medium. For adapters, self-locking.	Handgrip	13-311431	V

IMPRESSION AND LABORATORY ACCESSORIES

	Description	Code	REF	Price cat.
	Lab analog, with cone	IA STI	13-420347	B
	Transfer post, engages on the outer cone, 9 mm high Tighten with HT 1.25	TS STI	13-420345	B
	Transfer post, engages on the outer cone, 13 mm high Tighten with HT 1.25	TSX STI	13-420411	B
	Hex tool. Length 21 mm, Ø 1.25 mm	HT 1.25	13-425100	C
	Hex tool extra long. Length 45 mm, Ø 1.25 mm	HTX 1.25	13-425102	C
	Torx tool. Length 21 mm, Ø 1.25 mm	TT 1.25	13-425105	C
	Synthetic base castable, 7,4 mm high. For TLA STI and TCA STI	PS	13-420354	A
	Synthetic base castable, 4 mm high. For TLA15 STI	PS 15	13-420355	A
	Synthetic base castable, 4 mm high. For TLA20 STI	PS 20	13-420356	A
	Synthetic base castable, 4 mm high. For TLA25 STI	PS 25	13-420357	A

TPG® AND TPG® UNO SMALL TRAY



Size of closed tray
W 110 mm **L** 110 mm **H** 60 mm
 for all autoclaves

Autoclaveable up to 134 °C, not suitable for dry heat sterilizers.
 Contains all drills and tools for first works with the **TPG®** system.
 Material: autoclaveable plastic

Description

Insertion tool
 Insertion tool, length 20 mm
 Hex instrument, long
 Pathfinder drill
 Pathfinder drill
 Heatless drill for implants with conical core
 Heatless drill for implants with conical core
 Heatless drill for implants with conical core
 Countersink heads
 Standardized probe
 Twist drill
 Insertion tool, short 12 mm
 Insertion tool, medium 16 mm
 Insertion tool, long 20 mm
 Insertion tool for contra-angle, 30 mm
 Ratchet for all Hex tools and insertion tools

For system

TPG® & TPG® uno
TPG® & TPG® uno
TPG®
TPG® & TPG® uno
TPG® & TPG® uno
TPG® & TPG® uno
TPG®
TPG®
TPG®
TPG®
TPG®
TPG®
TPG®
TPG®
TPG®
TPG® & TPG® uno

Code

IT2 BCS
 IT K
 HT 1.25
 BCD 1
 BCDX 1
 DOS 2
 DOS 4
 DOS 5
 DRA 4.1 CS
 PDG
 Twist Drill 2.0
 IT1 TPG
 IT2 TPG
 IT3 TPG
 ITM TPG
 RAT 2

REF

13-900030
 13-462320
 13-425100
 13-900240
 13-900243
 13-455312
 13-455314
 13-455315
 13-425047
 13-425400
 13-90022
 13-425116
 13-425117
 13-425118
 13-425119
 13-425051

Price €






Small tray with TPG® content		13-60047-K	upon request
Small tray with TPG® uno content		13-60047-K	upon request
optional Adapter TPG for handgrip		13-425106	D
optional Adapter TPG for abutments OF5 STI, OF 5L STI		13-310530	D
optional Handgrip for adapters, self-locking		13-311431	V
Small tray with adapter and handgrip		13-60047-K	upon request



TPG® uno implants provide a sharp and cutting crestal thread, which engages in to the corticals of the extraction sockets. The apical compression thread provides perfect stability both in compressed spongy and in the 2nd cortical. **TPG® uno** implants may be used both in extraction sockets and in healed bone areas. They are used for circular bridges and segments. Under adequate loading conditions and if an excellent primary stability is achieved, the treatment can be performed in an immediate load protocol.

Due to their polished surface, **TPG®** and **TPG® uno** implants are extremely resistant against bacterial colonialization and peri-implantitis. **TPG®** and **TPG® uno** implants are not designed for single implant reconstruction.

TPG® uno IMPLANTS II POLISHED SURFACE, MATERIAL Ti6Al4V ELI, TRIPLE THREAD

		Description	Max. Ø	Max. endoss. Ø mm	Endosseous length	REF	Price cat.
TPG Uno 3.5 15		TPG Uno 3.5 12	4.25 mm	3.5	12	13-420120	K
		TPG Uno 3.5 15	4.25 mm	3.5	15	13-420121	K
		TPG Uno 3.5 17	4.25 mm	3.5	17	13-420122	K
		TPG Uno 3.5 19	4.25 mm	3.5	19	13-420123	K
TPG Uno 3.5 17		TPG Uno 3.5 21	4.25 mm	3.5	21	13-420124	K
		TPG Uno 3.5 23	4.25 mm	3.5	23	13-420125	K
		TPG Uno 3.5 26	4.25 mm	3.5	26	13-420126	K
		TPG Uno 3.5 29	4.25 mm	3.5	29	13-420127	K
TPG Uno 4.1 15		TPG Uno 4.1 12	4.95 mm	4.1	12	13-420130	K
		TPG Uno 4.1 15	4.95 mm	4.1	15	13-420131	K
		TPG Uno 4.1 17	4.95 mm	4.1	17	13-420132	K
		TPG Uno 4.1 19	4.95 mm	4.1	19	13-420133	K
TPG Uno 4.1 17		TPG Uno 4.1 21	4.95 mm	4.1	21	13-420134	K
		TPG Uno 4.1 23	4.95 mm	4.1	23	13-420135	K
<hr style="border-top: 1px dashed #000;"/>							
TPG Uno 3.5 15 SH		TPG Uno 3.5 12 SH	4.95 mm	3.5	12	13-420170	K
		TPG Uno 3.5 15 SH	4.95 mm	3.5	15	13-420171	K
		TPG Uno 3.5 17 SH	4.95 mm	3.5	17	13-420172	K
		TPG Uno 3.5 19 SH	4.95 mm	3.5	19	13-420173	K
		TPG Uno 3.5 21 SH	4.95 mm	3.5	21	13-420174	K
		TPG Uno 3.5 23 SH	4.95 mm	3.5	23	13-420175	K
		TPG Uno 3.5 26 SH	4.95 mm	3.5	26	13-420176	K
		TPG Uno 3.5 29 SH	4.95 mm	3.5	29	13-420177	K

Dimensions for all TPG uno

Neck height	3.8 mm
Neck diameter	2 mm (bendable up to 15°)
Abutment	3.9 w x 7.2 h (mm)

TPG® uno implants are delivered **incl. lab-set REF 13-462353** consisting of



Double analogue,
plastic

IA4/IAU

13-462111



Impression post castable,
internally edged

PA X

13-462136









Impression post castable,
internally round






TSPA 4

13-462029

INSERTION TOOL

	Description	Code	REF	Price cat.
 For single-piece TPG® uno	Adapter, for handgrips REF13-311430 & REF13-311431	Adapter Wst	13-310530	D
	Insertion tool short, for large head. Use with RAT2 and TW2. Length 7 mm	IT2 S BCS	13-900038	E
	Insertion tool medium, for large head. Use with RAT2 and TW2. Length 19 mm	IT2 BCS	13-900030	E
	Insertion tool for large head. Use with angle piece. Length 23 mm	IT2W	13-900039	C
	Torque wrench 10 - 70 Ncm. It is recommended to have the torque ratchets recalibrated by us once a year.	TW2	13-425402	S
	Ratchet for all hex instruments and insertion tools	RAT2	13-425051	K






IMPRESSION TAKING AND LABORATORY ACCESSORIES

	Titanium cap, radio opaque Ti6Al4V, weldable	MA5	13-462093	B
	Impression post castable, POM Internally round. Pack of 5	TSPA 5	13-462030	B
	Castable abutment Pack of 5	POB	13-462086	B
	Double analog, plastic Pack of 5	IA4/IAU	13-462111	B
	Double analog, metal	IA4/IAU	462112	A




HEATLESS® DRILLS DOS FOR IMPLANTS WITH CONICAL CORE


Surgical steel, color-coded, depth-coded and autoclaveable. The drill is marked with laser depth markings.
Use between 3,000 and 5,000 rpm with good cooling and intermittent drill technique.
Due to the extremely high cutting performance, you can work without pressure.



	Description	Colour	Max. working length	Code	REF	Price cat.
	DOS 1	yellow	17 mm	DOS 1	13-455311	D
	DOS 2	black	17 mm	DOS 2	13-455312	D
	DOS 3	red	17 mm	DOS 3	13-455313	D
	BCDX 1	yellow	15 mm	BCDX 1	13-900243	C
	BCDX 2	black	15 mm	BCDX 2	13-900244	C

DRILLS FOR THE PREPARATION OF THE 1ST CORTICAL

	Description		Code	REF	Price cat.
	C-Drill KM1 3.0 - 3.2	Cortical milling	C-Drill KM1	13-455300	E
	C-Drill KM2 3.7 - 4.1	Cortical milling	C-Drill KM2	13-455301	E
	C-Drill KM3 5.0	Cortical milling	C-Drill KM3	13-455302	E



We are certified DIN EN ISO 13485, and annex II of EEC Directive 93/42 EWG (2007).

The relative dimensions of parts on the pictures in this prospectus may differ from the reality (for technical reasons, or because items may have been further developed).

TPG® is a registered trade-mark. If implants are cleaned and re-sterilized in the dental office, the danger of creating infections is immanent, because in a dental office no processes for validated cleaning are available.

(The products of this catalogue are CE marked (class I) and CE 1936 marked (class IIa and IIb) according to 93/42/EC Directive).

Commercial products that are not monitored by our notified body are declared as third-party products.

Compilation and explanation of symbols on the packaging:



Batch No.



Sterilized by radiation



Catalogue number



Nonsterile



Intended for use by dentists or surgeons only



Single use product



Instruction for use



Expiry date



Store in a dry place



Store tightly keep closed



Do not use if packing is damaged



Do not re-sterilize



Manufacturer



Production date

Simpladent GmbH
Dorfplatz 11
CH - 8737 Gommiswald / SG
SWITZERLAND

Phone +41 (0)55 293 23 70
implants@simpladent.ch
www.simpladent-implants.com

PLEASE ASK YOUR LOCAL DISTRIBUTOR FOR THE VALID PRICE LIST

SIMPLADENT®



Dr. Ihde Dental GmbH
Erfurter Str. 19
D - 85386 Eching / München
Tel +49 (0)89 319 761 0
info@ihde-dental.de